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A wellness device for enhancing our own protective forces and correcting the natural electromagnetic radiation in natural bio-systems

DeVita AP
Model Mini

Instructions for use

Introduction

A unique device which has no analogues in the world. Each person needs to maintain a high rate of health in the modern world. The device is used to relieve the stressful burden caused by the presence of parasitic forms in the human body: viruses, bacteria, fungi, protozoa, helminths. And also the presence of endo- and exo-toxins that block the work of the main excretory systems, which is one of the main causes of reducing the state of health.

Method of work

Its action is based on the results of large-scale research by the scientist and engineer R. Rife and his supporters. Scientists have proved that every living form has its own electromagnetic radiation. Waves of a certain frequency, reproduced by the device, resonate with the own radiation of the microorganisms, enhancing their natural vibrations. Intensive fluctuations become disastrous for parasites. The complex AP Mini device has separate vibration frequencies for certain types of viruses, bacteria, helminths, protozoa and fungi, and, in addition, allows the deep purification of the body from toxins.

Features of the device

The device includes 14 basic complexes consisting of manual and automatic programs. Each user can create their own unique complexes (the whole device can hold up to 50) and has the ability to download new complexes and programs that will be created by the company's specialists.

Time of use

The DeVita AP Mini device should be used in the afternoon or at night, placing it under the pillow with the backside facing yourself.

To maintain a high wellness factor, to improve the state of health and for best adaptation, we recommend using the DeVita Energy Mini device in the morning and the DeVita Ritm Mini device during the day. The method of using DeVita devices is described in the complexes.

Complexes

- | | | |
|-------------------------------|--------------------------|--------------------|
| 1. General list of programs | 6. Active immune system | 11. Stop allergy |
| 2. Heart and circulation | 7. Movement without pain | 12. Men's health |
| 3. Digestion without problems | 8. Sport performers | 13. Women's health |
| 4. Stop stress | 9. Clean skin | 14. Body figure |
| 5. Stop insomnia | 10. Stop parasites | |

1. Complex «General list of programs»

- 1. Anti-Cold
- 2. Anti-Flu
- 3. Antivirus
- 4. Digestive tract health
- 5. Drainage antitox
- 6. Drainage, general
- 7. Help against poisoning
- 8. Intimate health
- 9. Mexican detoxification
- 10. No allergies
- 11. No bacteria
- 12. No candida
- 13. No chlamydia
- 14. No fungi
- 15. No flu, general
- 16. No helicobacter
- 17. No hepatitis A
- 18. Healthy lymph
- 19. No herpes
- 20. Without toxins
- 21. No herpes simplex
- 22. No herpes zoster
- 23. No mould & fungi
- 24. No nematodes & ascarides
- 25. No parasites, Digestive tract
- 26. No parasites, general
- 27. No protozoans
- 28. No schistosomes (blood-flukes) & trematods
- 29. No staphylococcus
- 30. No streptococcus
- 31. Weight Control
- 32. Mode 1 -automatic
- 33. Mode 2 -automatic
- 34. Mode 3 -automatic
- 35. Mode 4 -automatic
- 36. Cleansing -automatic
- 37. No parasites-automatic
- 38. No flu -automatic
- 39. Stop herpes S-automatic
- 40. Stop herpes Z-automatic
- 41. Without cold -automatic
- 42. Support immunity -automatic
- 43. Without helminths -automatic
- 44. Weight Reduction-automatic
- 45. Reduce Points -automatic

Description of automatic programs

Name of the program	Aim of program	Time
Mode 1 - automatic	Oppresses the activity of viruses, inhibits the vital activity of protozoal microorganisms. Optimal prevention of diseases when changing climatic zones, like visiting Asian countries.	6:47
Mode 2 - automatic	Removes toxic load caused by strepto- and staphylococcus. It is effective both for colds and for the prevention of coccal infections.	7:53
Mode 3 - automatic	Promotes the removal of bacterial complications of the gastrointestinal tract and urogenital system.	8:26
Mode 4 - automatic	It displays mycotoxins formed by molds and yeast fungi. Reduces the risk of accumulation of carcinogenic factors.	7:47

Cleansing - automatic	It stimulates the removal of toxins and the normalization of lymphatic drainage, strengthening the lymphatic vessels. Helps removing the excess fluid from the intercellular space, it has an antiseptic effect.	7:10
No parasites - automatic	Reduces the parasitic burden of the liver, pancreas, intestines, and also the toxic load caused by the vital activity of these forms..	7:38
No flu - automatic	It suppresses the activity of various strains of the flu virus, helps to relieve symptoms of intoxication and alleviate a person's condition during the disease..	6:38
Stop herpes S - automatic	Prevents the development of herpes simplex and its toxic effects on the body: the formation of vesicles on the mucous membranes, the skin and reduces itching.	7:35
Stop herpes Z - automatic	Reduces the activity of herpes zoster viruses, reduces the toxic signs of its manifestation on the skin surface and pain sensations.	5:47
Without cold - automatic	Prevents the development of respiratory viruses, reduces the symptoms of intoxication, leads to a decrease in inflammatory reactions and the likelihood of complications	8:09
Support immunity - automatic	Blocks the activity of a variety of viruses, their distribution in the human body, contributes to the activation of the immune system. Effective in preventing from infections during hypothermia and during the period of mass diseases.	5:03
Without helminths - automatic	Helps to impair the vital activity of helminths and their natural withdrawal from the body, reduces the toxic load. Effective in preventing from helminthiosis in children and adults.	3:18
Weight Reduction - automatic	Reduces the parasitic load and destroys viruses that cause excess weight. Normalizes the work of the autonomic nervous system, maintains a healthy tone of the body and improves the psychoemotional state.	4:28
Reduce Points - automatic	Clears the body of toxins, normalizes the lymph drainage, maintains immunity and increases vital energy.	4:31

Attention! If during the use of the programs of the first month nausea, headache, perspiration, or other symptoms of intoxication are observed, then additional application of this or other mode has to be carried out within the next 3 months. Each of them is involved in the course taking into account the life cycle of the infectious organism.

So:

- Modes 1-3: apply 14 days;
- Mode 4: 30 days.

The device is used daily at night

2. Complex «Heart and circulation»

Applied for problems with the cardiovascular and lymphatic systems.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month:

- 1st -2nd week — program «No schistosomes (blood-flukes) & trematods» daily use in the evening after the programs of DeVita Ritm Mini.
- 3rd -4th week — program «No chlamydia» daily use in the evening after the programs of DeVita Ritm Mini.
- 1-4 week — «Cleansing» daily at night.

3rd month: necessary use of the programs of the first month

4th month: «No parasites-automatic» daily at night time for 1 week.

3. Complex «Digestion without problems»

For chronic problems with the organs of the gastrointestinal tract and as an additional tool in case of poisoning.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month:

- 1st -2nd week — program «No schistosomes (blood-flukes) & trematods» daily use in the evening after the programs of DeVita Ritm Mini.
- 3rd -4th week — program «Drainage antitox» daily use in the evening after the programs of DeVita Ritm Mini.
- 1-4 week — «Cleansing» daily at night.

If necessary:

«No helicobacter» — daily in the evening for a month.

3rd month: necessary use of the programs of the first month

4th month: «No parasites-automatic» daily at night time for 1 week.

In case of poisoning:

- «Mode 3 - automatic»
- program «Help against poisoning»
Daily for 10 days, alternating with each other.

4. Complex «Stop stress»

In case of stress, fatigue, reduced performance, depressive states.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1- automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month: «Cleansing» daily at night.

3rd month: necessary use of the programs of the first month

4th month: «No parasites-automatic» daily at night time for 1 week.

5. Complex «Stop insomnia»

In case of a state of nervous excitement, trouble falling asleep, sleep disorders.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month: «Cleansing» daily at night.

3rd month: necessary use of the programs of the first month

4th month: «No parasites-automatic» daily at night time for 1 week.

6. Complex «Active immune system»

For the prevention of colds, as well as relief from bacterial and viral infections.

In the season of colds, in case of hypothermia:

- «Support immunity –automatic» 1-2 times a week

At the first symptoms of the common cold diseases (other than flu) during the first two days:

- «Without cold – automatic»

Then within 5 to 7 days alternating with each other:

- «Mode 2 - automatic»
- program «Anti-Cold»

At the first signs(symptoms) of the flu, for the first 2 days alternating with each other:

- program «Antivirus»
- «No flu –automatic»

Further, daily for 5 days: «No flu – automatic»

Upon completion, use «Mode 2 - automatic» daily for 5-7 days, also if necessary the program «No bacteria».

7. Complex «Movement without pain»

In case of problems with the musculoskeletal system: pain, inflammation of the joints, muscles, spine or injuries.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month:

- «Cleansing» daily at night.

- program «No bacteria» and «Drainage antitox». Daily in the evening after the DeVita Ritm Mini programs for 14 days.

3rd month: necessary use of the programs of the first month

4th month: «No parasites-automatic» daily at night time for 1 week.

8. Complex «Sport performers»

With intense physical activity, injuries, as well as to reduce the toxic load and a fast recovery after workouts.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month: «Cleansing» daily at night.

3rd month: necessary use of the programs of the first month

4th month: «No parasites-automatic» daily at night time for 1 week

9. Complex «Clean skin»

With inflammatory processes on the skin, fungal burdening of hair and nails, herpetic weights (shingles).

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month: «Cleansing» daily at night.

If necessary:

- in case of fungal burdening program «No fungi» daily in the evening after the programs of DeVita Ritm Mini

3rd month:

- necessary use of the programs of the first month
- in case of herpetic burdening: «Stop herpes S-automatic» or «Stop herpes Z-automatic» at night time.

4th month: «No parasites-automatic» daily at night time for 1 week.

10. Complex «Stop parasites»

To reduce health risks of: pet owners, gardeners, travelers, families with children under 12 years. In case of appetite changes, unexplained pain, inflammation, appearance of allergies.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month: «Cleansing» daily at night.

3rd month: necessary use of the programs of the first month

4th month: «No parasites-automatic» and «Without helminths – automatic», alternating with each other every day at night for 1 month.

Additionally:

In case of cough, program «No bacteria» — daily in the evening after the programs of DeVita Ritm Mini, for 14 days.

In case of herpetic weights:

- program «No herpes simplex» at the first sign of a rash every 3 hours, then every day 1-2 times a day «Stop herpes S-automatic» at night for 1 month.
- program «No herpes zoster» at the first sign of a rash every 3 hours, then every day 1-2 times a day «Stop herpes Z-automatic» at night for 1 month.

When visiting countries in Asia and Africa, to prevent viral load of the liver, use the program — «No hepatitis A» — daily 1-2 times a day for a month.

11. Complex «Stop allergy»

With various allergic reactions to: food, dust, pets, pollen of plants.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

Additionally: program «No allergies» in the evening after the programs of DeVita Ritm Mini.

2nd month:

- 1st -2nd week — program «No schistosomes (blood-flukes) & trematods» daily use in the evening after the programs of DeVita Ritm Mini.
- 3rd -4th week — program «No allergies» daily use in the evening after the programs of DeVita Ritm Mini.
- 1-4 week — «Cleansing» daily at night.

3rd month:

- necessary use of the programs of the first month
- in case of herpetic burdening: «No allergies» daily at night.

4th month:

- «No parasites-automatic» daily at night time for 1 week.
- program «No allergies» in the evening after the programs of DeVita Ritm Mini.

Attention! In case of pollinosis, the procedure should be started 2 months before the flowering of plants.

A stable positive result for solving problems with allergies can be obtained with regular use of the programs of the complex for 6 months.

12. Complex «Men's health»

In case of psychological and functional problems with the organs of the male urogenital system.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month:

- «Cleansing» daily at night.
- program «No chlamydia» in the evening after the programs of DeVita Ritm Mini.

3rd month: necessary use of the programs of the first month

4th month: «No parasites-automatic» daily at night time for 1 week

13. Complex «Women's health»

In case of functional disorders of the female urogenital system, for harmonizing the work of the organism in premenstrual and menopausal period.

1st month of daily use, during night time, according to the following scheme:

- 1st week — «Mode 1 - automatic»
- 2nd week — «Mode 2 - automatic»
- 3rd week — «Mode 3 - automatic»
- 4th week — «Mode 4 - automatic»

2nd month:

- «Cleansing» daily at night.
- program «No chlamydia» in the evening after the programs of DeVita Ritm Mini.

3rd month: necessary use of the programs of the first month

4th month: «No parasites-automatic» daily at night time for 1 week

14. Complex «Body figure»

For soft and comfortable weight loss, regulation of appetite.

Attention! To obtain the most effective result before using this complex, we recommend first to reduce the possible parasitic and toxic load, using the basic programs («Mode 1-4 automatic», «Cleansing», «No parasites-automatic») of the complex «Stop parasites» within 3-4 months.

4th-6th month: programs «Weight Reduction-automatic» and «Reduce Points – automatic» every day, alternating with each other at night

For detailed information on the complexes, you can visit the company's website www.deholding.info under the section «Products» – «Complexes».

With regular, correct use of DeVita Mini devices any user will receive a unidirectional positive result of effective recovery, youth and longevity.

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Deta-Elis Europa GmbH